

Fiber-optic pH Minisensors



For fast pH measurements in small sample volumes


FireSting[®]-PRO Read-out Meters

with 1-, 2- or 4-channels for connecting sensors



Minisensors for different pH ranges

- pH 7.0 - 9.0 (-PK8)
- pH 5.5 - 7.5 (-PK6.5)
- cable length of ca. 2 m
- tip diameter of 500 μm
- fixed-needle type or bare fiber
- response time t_{90} at 25°C: <10 s
- for screening, single-use and fast measurements in small volumes



Ideal for fast pH measurements in small samples!